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Bartosz Machoński

**MalTec Health Shoulder Massager
PRO-Comfy30W-02
Instruction Manual**

Before usage please read whole instruction carefully.

Warning:

To reduce causing a fire, burns, electric shock or any other injury please observe the following instructions.

- Once finished using the device please pull out the plug from the socket.
- To avoid overheating the product please do not use under pillows, blankets etc...
- Please keep away from children.
- Please keep away from water and liquids.
- Users with physical impairment or a disability should use under guidance and after reading the instructions.
- Please use this product according to this manual provided.

Danger:

In order to reduce the risk of electric shock, please follow the instructions below:

- After using or before cleaning product, please pull out the plug.
- If the product falls into the water, don't touch the wet product. You must pull out the plug.
- Do not use the product when taking a shower or a bath.
- Do not place the product in water or other liquid.

Troubleshooting

If you encounter problems when using this device please refer to the following tips and solutions. If the massager is still in a faulty state or not available, please turn off the power and unplug the device. User should not disassemble or repair the massager themselves, please revert back to the seller for guidance on repair or replacement.

Condition	Possible reasons/Solutions
Unable to start the massager	Make sure the power plug is in socket. Confirm the switch is on. Confirm if the button has been pressed.
When using, it suddenly stops to work.	Please check whether the power plug is pulled out or not. Confirm if the plug has a bad connection. Automatic timer sets for 15 minutes, time has passed. If the overheating protection device starts, please turn off the main power supply, and wait for the massage cooling. If the problem persists, please do not continue to use and contact the company repair service centre.

Regular massage can improve the whole body function.

Improve the quality of sleep.

Massage can make the body relax, eliminate fatigue, and relieve pressure. It also can improve sleep, and improve sleep quality

Improve body function

Massage can promote blood circulation for the body to provide adequate nutrients, improve lymph detoxification.

Condition gastrointestinal function

Massage can improve gastrointestinal function, promote the new supersedes the old, improve digestion.

Help digestion, relieve constipation and other gastrointestinal adverse circumstances.

Massage Neck & Shoulder

Exquisite curve design allows a tight massage for the neck and shoulder.

Kneading massager

This is a deep massage. The massage heads knead and massage upward or downward.

Warming Function

Warming function this ease muscle tension and pain.

Home usage or on the go

Equipped with an external adapter, you can also use it by 12V vehicle power supply.

Notice

- The transformer (DC12V 2A) or rated voltage must be provided.
- Ensure the device and the transfer work together.
- Don't remove the connection by force this can cause damage.
- Device should not be placed on a hard surface, this can cause damage to the massage and heating functions.

Disposal of electrical equipment

The symbol of the crossed-out wheeled waste container placed on the products informs that used electrical and electronic equipment (hereinafter also "waste equipment" or "electro-waste") may not be disposed of with municipal waste. By throwing e-waste into the container, you pose a threat to the environment. Used equipment may contain substances (e.g. lead, cadmium, chromium, bromine, mercury, freon) that are dangerous to human health and life and to the environment. By sorting and handing over used equipment for processing, recovery, recycling and disposal, you protect the environment against pollution and contamination, you contribute to reducing the use of natural resources and lowering the production costs of new equipment.

CE declaration

Malis B. Machoński sp. k. Declares that the MalTec Health PRO-Comfy30W-02 shoulder massager meets all the requirements of the EC 2014/30/EU directive and meets the following standards:

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