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PL

MASAŻER DO RĄK
HANDPULSE BLACK

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INSTRUCTION MANUAL EN

HAND MASSAGER
HANDPULSE BLACK

**PLEASE READ THIS USER MANUAL BEFORE USING THE DEVICE!
FOLLOW ALL SAFETY INSTRUCTIONS TO AVOID DAMAGE CAU-
SED BY IMPROPER USE! KEEP THE USER MANUAL FOR FUTURE
REFERENCE.**

SAFETY INSTRUCTIONS

Please follow the doctor's advice for the following patients:

The patients who have had hand surgery, swelling and inflammation, congestion, severe arthritis, trauma and other hand diseases;
People with heart, brain and lung diseases;
Blood diseases with bleeding tendency;
Osteoporosis and soft tissue diseases, stroke;
Infants and pregnant women.

1. Use Type-C plug connection.(Note the shape of the insertion port)
2. Before use, Please make sure that you have no jewelry(watch ring, Bracelet ring, etc.) in your hands, so as to avoid sharp objects piercing the airbag or squeezing the airbag and causing jewelry damage.
3. It is forbidden to use for a long time. Normally, massage time should not exceed 30 minutes.
4. Please do not use this product while driving.
5. Please do not use it after getting drunk or taking sleeping pills.
6. Please use DC5V1~2.5A charger.
7. This product is specially designed as household electronic health care products. It is not suitable for diagnosis, treatment and other professional or commercial purposes.
8. Please keep this product away from water, fire and corrosive environment.
9. Please do not put sharp objects into this product.
- 10 Please do not use when hands are wet.
11. If the power cord is damaged, stop using it immediately. Contact the product service center for maintenance. Do not disassemble the product at will.
12. If you feel unwell, please stop using it in time.
13. Please do not wrap products with towels, quilts, clothes, etc. when you use it.

PRODUCT DESCRIPTION

USER INSTRUCTIONS

1. Long press the power key to turn on the machine, default mode 1 (Gentle mode), strength 1 gear, the corresponding mode indicator light strength indicator light are on. Long press the power key again to turn off the machine.
2. Touch the heating key, the corresponding heating indicator lights up, and the heating function is turned on at level 1. Touch it again, the corresponding heating lights up, and the heating function is turned on at level 2. Touch it again, the corresponding heating indicator lights up, and the heating function is turned on at level 3. Touch it again, the heating indicator goes out, and the heating function is turned off.
3. Press the vibration key, the corresponding vibration indicator lights up, and the vibration function is turned on at level 1. Press once more, the corresponding vibration lights up, and the vibration function is turned on at level 2. Press once more, the corresponding vibration indicator lights up, and the vibration function is turned on at level 3. Press once more, the vibration indicator goes out, and the vibration function is turned off.
4. Press the mode key, can adjust the mode, press once, the corresponding mode indicator lightlights up, soothing mode. Press again, the corresponding mode indicator light up, powerful mode.
5. Press the pressure key to adjust the pressure, there are 5 gears, 1-2-3-4-5, the corresponding indicator lights are on
6. Press the timer key to adjust the timer time, 10-15-20 minutes. The countdown begins, time is up, the machine shut down.
7. When charging, the number will display 90 percent of the battery, and when fully charged, it will display the number 100.
8. When the power is on, the buzzer will beep three times and the indicator light will flash several times, indicating that the power is exhausted and the product needs to be charged.

- Do not use while charging.
- When the product is charging, the power button flashes until it is fully charged.
- After applying the hand cream, massage with plastic gloves can help absorb hand cream and repair skin.
- Please don't use products with heating function for a long time to prevent low temperature burns.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved

MAINTENANCE AND CLEANING / STORAGE

1. Use towel to remove dust and other dirt, and place it in a dry and ventilated place.
2. Do not place in high temperature, humidity Maintenance Soft tissue diseases; stroke; infants and pregnant women; and storage and direct sunlight.
3. When the storage period exceeds 6 months, it should be checked out and power-on.
4. When plugged into the USB socket, the machine can be used. But we strongly recommend against it.

TROUBLESHOOTING

If there is a malfunction during the use of this product, please check it as follows. If the fault can not be eliminated or the product still cannot run, please turn off the power and contact the product service center for maintenance.

Basic Faults	Possible causes	Troubleshooting Measures
The product can't start working	Check the remaining power of the battery The power may be weak.	Charging
The product cannot be re-charged	Is the power adapter connected to the data line correctly?	Re-plug until LCD power display flickers
	Is the data line and USB socket intact?	Check the socket If there is any damage, please contact the manufacturer
	Are the specifications of adapters and data lines in line with operational requirements?	View specifications and proofread the technical parameters provided in the instructions
The product can be turned on, but can't work	The remaining power may be weak.	Charging

TECHNICAL DATA

RATED VOLTAGE	DC 5V
BATTERY CAPACITY	2500mAh
RATED POWER	6W
WORKING TIME	10/15/20 min
CHARGING	5V 1A
WEIGHT	0.85kg
SIZE	267*162*120mm

WARRANTY

The seller is liable for defective goods for 2 years from the date of delivery.

All pictures and drawings are for reference only.

Disposal of electrical equipment

The symbol of the crossed-out wheeled waste container placed on the products informs that used electrical and electronic equipment (hereinafter also "waste equipment" or "electro-waste") may not be disposed of with municipal waste. By throwing e-waste into the container, you pose a threat to the environment. Used equipment may contain substances (e.g. lead, cadmium, chromium, bromine, mercury, freon) that are dangerous to human health and life and to the environment. By sorting and handing over used equipment for processing, recovery, recycling and disposal, you protect the environment against pollution and contamination, you contribute to reducing the use of natural resources and lowering the production costs of new equipment.

DECLARATION OF CONFORMITY

Malis B. Machoński Sp. K. declares that device described in this manual meets all the requirements of the:

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