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Bartosz Machoński

**MalTec Health
Pillow Massager
PRO-Comfy30W-01
Instruction Manual**

Thanks for you to select our product. Before use, please carefully read the manual and pay more attention to safety notice items to correctly operate and use.

Please take good care of the manual.

Safety Items

Use Condition

- Do not use in wet or dirty environment to prevent any trouble.
- Do not use if room temperature is above 104 degrees F (40°C).
- Do not use the product outdoor.

Precautions

If you have the following disorder or condition, please consult your physician before use: Heart Failure, Kidney Failure, Infection of the superficial veins (called phlebitis) or soft tissue (called cellulitis) in the legs or elsewhere, Blood clots in the legs and Bleeding disorders.

Unfit people to use this product. If you have the following state, before use, please consult the doctor.

- Curer who is receiving the cure from the doctor or who is feeling bad.
- Malignancy Patient.
- Pregnant Woman or Woman in menses.
- Urgent illness Patient.

Notice

During massaging, the massage head should be near two sides of the neck. Otherwise, it might result in an ache.

Operating Safety Measures

- Do not operate with wet hands.
- Do not operate in water.
- Do not drag the controller cable.
- Do not place heavy objects on the controller.

Massager Operations

- During massage, ensure that the massage heads are on equal sides of the neck to eliminate unnecessary discomfort.
- Do not sleep while using the massager.
- Do not use while intoxicated.
- Do not directly massage hand, elbow parts, knee, breast or abdominal area.
- Do not exceed 30 minutes of use.

Use Method

1. Plug into power socket or use car power adaptor while in a vehicle.
2. Turn on switch key, massage pillow enter into massaging state. Advisable to use in 10 to 15 minutes increments.
3. Stop continuing to use, please turn off the switch button. Then plug out from power socket or power adapter

Frequently Asked Questions

Sequence No	Phenomenon	Possible Reason	Solution
1	In operating app hear „Xiu xiu dada” sound. Product in working order	Or the frictional sound between massage ball and fabric packing	No need any disposal
2	Open to fail	Product is too hot to protect	Turn off for 30 minutes to use it again
		Power Plug or Smoke source socket is not plugged well	Plug the socket well
		Direction Switch dial to middle position	Direction toward left or right

If dispose the above, still there is any abnormal situation. Please contact with local distributor and request them to maintain and repair.

Stop using

When there is the following, please immediately stop using the product. Close power Switch and Plug out of the plug.

- When suddenly splash water to this product.
- The relevant parts will be broken or inner parts will expose.
- During massaging find any abnormal state such as fume, peculiar smell, etc.
- Suddenly power broken.
- Thunder occur.

Safety Items

Product Assembly & Maintenance

- Please don't privately disassembly the packing of this product.
- Do not touch the inner parts of this product to avoid any trouble.
- This product can be maintained only by the professional service worker, do not privately disassemble or maintain.

About Power Adapter & Car Power Cable

- Please check whether the power adapter accord with the specified voltage, current, and matches the output.
- When plugging out power cable hold power plug, NOT power cable.
- Power cables don't be pushed down by other heavy goods.

Inner Part Name & Functions Appearance Construction

Disposal of electrical equipment

The symbol of the crossed-out wheeled waste container placed on the products informs that used electrical and electronic equipment (hereinafter also "waste equipment" or "electro-waste") may not be disposed of with municipal waste. By throwing e-waste into the container, you pose a threat to the environment. Used equipment may contain substances (e.g. lead, cadmium, chromium, bromine, mercury, freon) that are dangerous to human health and life and to the environment. By sorting and handing over used equipment for processing, recovery, recycling and disposal, you protect the environment against pollution and contamination, you contribute to reducing the use of natural resources and lowering the production costs of new equipment.

CE declaration

Malis B. Machoński sp. k. Declares that the MalTec Health Pillow Massager PRO-Comfy30W-01 meets all the requirements of the EC 2014/30/EU directive and meets the following standards:

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