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Bartosz Machoński

**MalTec Health  
Pillow Massager  
ULTRA Comfy 30W  
Instruction Manual**

Thanks for you to select our product. Please read the manual for safety precautions and enjoy your new massager. Keep this manual convenient and available for reference.

## **Safety Items**

- Do not use in wet or dirty environment.
- Do not use if room temperature is above 104 degrees F (40°C).

## **Precautions**

If you have the following disorder or condition, please consult your physician before use: Heart Failure, Kidney Failure, Infection of the superficial veins (called phlebitis) or soft tissue (called cellulitis) in the legs or elsewhere, Blood clots in the legs and Bleeding disorders.

## **Operating Safety Measures**

- Do not operate with wet hands.
- Do not operate in water.
- Do not drag the controller cable.
- Do not place heavy objects on the controller.

## **Massager Operations**

- During massage, ensure that the massage heads are on equal sides of the neck to eliminate unnecessary discomfort.
- Do not sleep while using the massager.
- Do not use while intoxicated.
- Do not directly massage hand, elbow parts, knee, breast or abdominal area.
- Do not exceed 30 minutes of use.

## **Points For Attention**

1. To prevent electrical shock, DO NOT pull manual control knob. Disconnect power cord by the plug. Do not pull at the plug body.
2. If massager does not function, turn power off and discontinue use.
3. Do not allow children to use or play with the massager.
4. Discontinue use in the presence of discomfort.
5. Seek medical advice prior to use if you are pregnant, have a fever, malignant tumour, have heart disease, recovering from a surgical procedure, or have an acute disorder.
6. Do not use in places with high temperature, high pressure, near water or fire sources.
7. Keep away from sharp items. Repair if product's material is cut or torn.
8. Disconnect from power source before wiping with damp cloth.
9. Do not use if power cord is frayed or damaged.

## To Operate

1. Plug into power socket or use car power adaptor while in a vehicle.
2. Dial direction switch key to start massage pillow. Advisable to use in 10 to 15 minutes increments.
3. Discontinue by turning dial to the mid position and unplugging power cord or car power adaptor.

## Frequently Asked Questions

Sequence No	Phenomenon	Possible Reason	Solution
1	In operating app hear „Xiuxiu dada” sound. Product in working order	Or the frictional sound between massage ball and fabric packing	No need any disposal
2	Open to fail	Product is too hot to protect	Turn off for 30 minutes to use it again
		Power Plug or Smoke source socket is not plugged well	Plug the socket well
		Direction Switch dial to middle position	Direction toward left or right

If dispose the above, still there is any abnormal situation. Please contact with local distributor and request them to maintain and repair.

## CE declaration

Malis B. Machoński sp. K. Declares that the device MalTec Health Massager ULTRAComfy30W pillow meets all the requirements of the Directive LVD 2014/35 / EU, EC 2014/30 / EU and meets the following standards:

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